









Littmann

Stethoscopes

Choose the right stethoscope to help you perform at your best

	CORE Digital	Master Cardiology	Cardiology IV	Classic III	Classic II Pediatric and Infant	Lightweight II S.E.
	<p>Switch from analog listening to digital. Connects to Eko software on a smart device** to capture, store and share data.</p> <p>Chestpiece diameter double-sided Adult: 1.7"/4.3cm Pediatric: 1.4"/3.4cm</p> 	<p>Best acoustic response in Littmann mechanical line.</p> <p>Chestpiece diameter single-sided: 1.9"/4.8cm</p> 	<p>Chestpiece diameter double-sided Adult: 1.7"/4.3cm Pediatric: 1.4"/3.4cm</p> 	<p>Chestpiece diameter double-sided Adult: 1.7"/4.3cm Pediatric: 1.4"/3.4cm</p> 	<p>Specially designed for the smallest patients.</p> <p>Chestpiece diameter double-sided Pediatric: 1.4"/3.4cm, 1.1"/2.8cm Infant: 1.1"/2.8cm, 0.7"/1.7cm</p> 	<p>Teardrop-shaped chestpiece fits under blood pressure cuff.</p> <p>Chestpiece diameter double-sided 2.0"x1.8"/5.2cm x 4.5cm 1.4"x1.3"/3.5cm x 3.2cm</p> 
Acoustic performance level	10+	10	9	8	7	6
Length	27"/68cm	27"/68cm	27"/68cm 22"/56cm	27"/68cm	28"/71cm	28"/71cm
Weight	200g+	161g-200g	161g-200g	131g-160g	<100g, 100g-130g	100g-130g
Warranty	1-3 years	7 years	7 years	5 years	3 years	2 years
Diaphragm tunable	●	●	●	●		●
Up to 40x amplification*	●					
Noise cancellation	●					
Mobile connectivity	●					
Dual-lumen tubing	●	●	●			
Additional special finishes	●	●	●	●	●	
Convertible open bell	●	●	●	●	●	●
Soft sealing eartips	●	●	●	●	●	●
Tunable, dual-sided chestpiece	●		●	●		●

*At peak frequency vs. analog mode. Updated November 2023.
 **Subscription may be required.

© Solventum 2024. Solventum, the S logo and Littmann are trademarks of Solventum or its affiliates. 70-2013-1757-8